




Model # 800 Stage

*Shown with Optional
Side load paddle
w/rotator #WSR-800*



PROPRIETARY
This document contains information that is the sole property of Franklin Mechanical & Control Inc. and cannot be copied, reproduced or utilized in any way without the express written consent of Franklin Mechanical & Control Inc.

- **9 lbs. 8 oz.**
- **Rack and Pinion Drive**
- **Smooth Ball Bearing Travel**
- **Custom Holes Available**

- **Black Anodized**
- **Stainless Steel Raceways**
- **One Year Warranty**
- **Co-Axial Knobs**

Franklin Mechanical & Control, Inc.
8901-B Murray Ave., Gilroy, CA 95020 408-842-0313 www.franklinstages.com