

**STEREO MICROSCOPE REFERENCE CHART**
[www.oem-optical.com](http://www.oem-optical.com)
**STANDARD Features / Options:**

Make	Model	Type	Zoom Ratio	Body Magnification Range	Maximum Magnification (w/aux lenses)	Standard Working Distance W.D	W.D (w/aux lenses)	STANDARD Features / Options:								
								Photo-port Available	Optional eyepieces	Optional objectives	Standard Eyetube Angle	Ergonomic Head Option	Modular Infinity Design	MFG Country	PDF	
<b>OEM-OPTICAL stereozoom</b>																
OEM	<a href="#">SZ445</a>	Stereozoom	4.3:1	.7x - 3x	3.5x - 120x	97mm	97 - 160mm	•	•	•	45°		China	<a href="#">SZ445</a>		
OEM	<a href="#">SZ645</a>	Stereozoom	6.3:1	.8x - 5x	4x - 200x	115mm	43 - 211mm	•	•	•	45°		China	<a href="#">SZ645</a>		
OEM	<a href="#">Classic SZ-4</a>	Stereozoom	4.3:1	.7x - 3x	3.5x-60x	100mm	100 - 168mm		•	•	30°		China	<a href="#">Zoom 4</a>		
OEM	<a href="#">SZ600</a>	Stereozoom	6.3:1	.8x - 5x	2.4x - 250x	70mm	32 - 124mm	•	•	•	20°	•	•	China	<a href="#">SPZ</a>	
OEM	<a href="#">SZ800</a>	Stereozoom	8:1	.8x - 6.5x	2.4x - 325x	70mm	32 - 124mm	•	•	•	20°	•	•	China	<a href="#">SPZ</a>	
OEM	<a href="#">SZ1000</a>	Stereozoom	10:1	.8x - 8x	2.4x - 400x	70mm	32 - 124mm	•	•	•	20°	•	•	China	<a href="#">SPZ</a>	
<b>Meiji Stereozoom</b>																
MEIJI	<a href="#">RZ</a>	Stereozoom	10:1	1x - 10x	5x - 400x	76mm	29 - 164mm	•	•	•	45°	•	•	Japan	<a href="#">RZ</a>	
MEIJI	<a href="#">EMZ-1</a>	Stereozoom	3:1	1x - 3x	3x - 180x	93mm	34 - 251mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-2</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	93mm	34 - 251mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-2TR</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	93mm	34 - 251mm	•	•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-5</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	93mm	34 - 251mm	•	•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-6</a>	Stereozoom	6.4:1	.62x - 4x	2.1x - 270x	108mm			•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-8U</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	104mm	38 - 300mm	•	•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-9</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	93mm	34 - 251mm		•	•	60°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-10</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	110mm	45 - 326mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-12</a>	Stereozoom	6.2:1	.4x - 2.5x	2.5x - 75x	185mm	185 - 287mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-12TR</a>	Stereozoom	6.2:1	.4x - 2.5x	2.5x - 75x	185mm	185 - 287mm	•	•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-13</a>	Stereozoom	7:1	1x - 7x	3x - 420x	90mm	32 - 230mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-13TR</a>	Stereozoom	7:1	1x - 7x	3x - 420x	90mm	32 - 230mm	•	•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">EMZ-13VX</a>	Stereozoom	7:1	1x - 7x	3x - 420x	90mm	32 - 230mm		•	•	45°			Japan	<a href="#">EM series</a>	
MEIJI	<a href="#">Z-7100</a>	Stereozoom	6.4:1	.7x - 4.5x	2.1x - 270x	110mm	45 - 332mm		•	•	0°			Japan	<a href="#">EM series</a>	
<b>Labomed Stereozoom</b>																
LABOMED	<a href="#">CZM-4</a>	Stereozoom	4:1	1x - 4x	5x - 160x	100mm	25 - 150mm	•	•	•	45°		India	<a href="#">CZM4</a>		
LABOMED	<a href="#">CZM-6</a>	Stereozoom	7:1	0.65x - 4.5x	3.5x - 180x	100mm	25 - 150mm	•	•	•	45°		India	<a href="#">CZM6</a>		
LABOMED	<a href="#">CSM-2</a>	2 step turret	3:1 or 2:1	1x/3x, 2x/4x	10x - 80x	128 , 95mm	77 - 128mm		•		45°		India	<a href="#">CSM-2</a>		