

Penguin/Pro-FS



Foot Switch for Hands Free Image Capture

- The Penguin/Pro-FS Foot Switch allows for image capture while focusing your attention on microscope operation
- A Great help for applications that require handsfree operation
- Works with both Penguin & Pro Series Camera Systems
- Connects via mini DIN directly to your Penguin PCI card
- Durable and rugged design for years of everyday use
- 10ft (3m) cable



Pixera Corporation (USA/International)

Tel: (1) 408 341 1800, Tollfree: 1-888-474-9372 Fax: (1) 408 341 1818 sales@pixera.com www.pixera.com

Pixera UK Ltd (Europe)

Tel: +44 (0) 118 931 9635 Fax: +44 (0) 118 931 9636 eusales@pixera.com

Pixera Japan (Japan)

Tel: +81 44 850 2800 Fax: +81 44 850 2868 jpsales@pixera.com

© 2002 Pixera Corporation. All rights reserved. All trademarks are the trademarks of their respective owners. Specifications are subject to change without further notice.

Technical Specification

Ergonomics:

3.45"(W)x3.98"(D)x1.28"(H, Uncompressed)

Weight: 14.7oz (0.4kg) (Without cable)

Cable: Standard 7 pin Mini DIN Connector:

Environmental: · Corrosion proof reinforced polymeric material for rugged duty

Non-skid basepad and raised ribs for positive traction during operation