

Penguin/Pro-CS

edoc

Camera Switch

- Operate 2 camera systems from ONE computer
- Works with all Penguin/Pro Series camera systems
- Great for users of both color and monochrome camera systems
- Save table-top space by using ONE computer for two camera systems
- Durable electro-mechanical design
- Small size for table-top use
- Works with both Windows® and Mac® operating systems
- Included: Penguin/Pro-CS switch unit and 10ft (3m) cable



Pixera Corporation (USA/International)

Tel: (1) 408 341 1800, Tollfree: 1-888-474-9372 Fax: (1) 408 341 1818 sales@pixera.com www.pixera.com

Pixera UK Ltd (Europe)

Tel: +44 (0) 118 931 9635 Fax: +44 (0) 118 931 9636 eusales@pixera.com

Pixera Japan (Japan)

Tel: +81 44 850 2800 Fax: +81 44 850 2868 jpsales@pixera.com

© 2002 Pixera Corporation.

All rights reserved.

All trademarks are the trademarks of their respective owners. Specifications are subject to change without further notice.

Technical Specification

 Type:
 Electro-mechanical switch (No separate power needed)

 Size:
 5.9" (W) x 3.9" (D) x 2.16" (H)

 150mm(W) x 100mm(D) x 55mm(H)

 Weight:
 1Lbs 2.6oz (0.53kg) (Without cable)

 Connectors (44-pin):
 1. Camera A

Camera B
 I/O Connector (PCI)

Cable: 10ft (3m)
Color: Off White, textured surface

System Layout

