

**MICROSCOPE LIGHTING PRODUCTS REFERENCE CHART**

[www.oem-optical.com](http://www.oem-optical.com)

STANDARD Features / Options:

Make	Model	Type	Approvals	Illumination	Output	Lamp	Input Power	STANDARD Features / Options:							
								Aperture Control	Remote Analog Control	Electro-Shutter Control	Filter Holder	5/8" Fiber Output	1" Fiber Output	MFG Country	PDF
<b>FIBER OPTICS LIGHT SOURCE:</b>															
Techniquip	<b>21AC Series</b>	Fiberoptic		Halogen	21v 150watt	<a href="#">EKE/EJV</a>								USA	<a href="#">21AC spec</a>
	<a href="#">21AC1K100X</a>	" "		" "	" "	" "	110 VAC			Opt				USA	
	<a href="#">21AC1K110X</a>	" "		" "	" "	" "	" "	•		Opt				USA	
	<a href="#">21AC1K200X</a>	" "		" "	" "	" "	" "			Opt	•			USA	
	<a href="#">21AC1K300X</a>	" "		" "	" "	" "	" "			Opt	•	•		USA	
	<a href="#">21AC1K310X</a>	" "		" "	" "	" "	" "	•		Opt	•	•		USA	
	<a href="#">21AC2K100X</a>	" "	CE	" "	" "	" "	220 VAC			Opt				USA	
	<a href="#">21AC2K110X</a>	" "	CE	" "	" "	" "	" "	•		Opt				USA	
	<a href="#">21AC2K200X</a>	" "	CE	" "	" "	" "	" "			Opt	•	•		USA	
	<a href="#">21AC2K300X</a>	" "	CE	" "	" "	" "	" "			Opt	•			USA	
	<a href="#">21AC2K310X</a>	" "	CE	" "	" "	" "	" "	•		Opt	•			USA	
Techniquip	<b>21DC Series</b>	Fiberoptic	CE/UL	Halogen	21v 150watt	<a href="#">EKE/EJV</a>								USA	<a href="#">21DC spec</a>
	<a href="#">21DCUK100X</a>	" "	CE/UL	" "	" "	" "	Auto 100-240vac		•	Opt				USA	
	<a href="#">21DCUK110X</a>	" "	CE/UL	" "	" "	" "	" "	•	•	Opt				USA	
	<a href="#">21DCUK200X</a>	" "	CE/UL	" "	" "	" "	" "		•	Opt		•		USA	
	<a href="#">21DCUK300X</a>	" "	CE/UL	" "	" "	" "	" "		•	Opt	•			USA	
	<a href="#">21DCUK310X</a>	" "	CE/UL	" "	" "	" "	" "	•	•	Opt	•			USA	
Techniquip	<b>FOI-150 Series</b>	Fiberoptic		Halogen	21v 150watt	<a href="#">EKE/EJA</a>								USA	<a href="#">FOI spec</a>
	<a href="#">FOI-150</a>	" "		" "	" "	" "	110 VAC			Opt	•			USA	
	<a href="#">FOI-150-220</a>	" "		" "	" "	" "	220VAC			Opt	•			USA	
	<a href="#">FOI-150-BREM</a>	" "		" "	" "	" "	110VAC	•		Opt	•			USA	
Techniquip	<b>FOI-250 Series</b>	Fiberoptic		Halogen	24v 250watt	<a href="#">ELC</a>								USA	<a href="#">FOI spec</a>
	<a href="#">FOI-250</a>	" "		" "	" "	" "	110 VAC			Opt	•			USA	
<b>FLUORESCENT RINGLIGHTS:</b>															
Techniquip	<a href="#">TLC-5</a>	Ringlight		Fluorescent	Standard Ballast	<a href="#">WTL</a>	110 or 220VAC							USA	<a href="#">TLC5 spec</a>
Techniquip	<a href="#">TLC-7</a>	Ringlight		Fluorescent	HF - Variable	<a href="#">XTL</a>	Auto 100-240vac							USA	<a href="#">TLC7/8</a>
Techniquip	<a href="#">TLC-8</a>	Ringlight		Fluorescent	HF - Fixed	<a href="#">XTL</a>	Auto 100-240vac							USA	<a href="#">TLC7/8</a>
<b>LED RINGLIGHTS:</b>															
Techniquip	<a href="#">AN-2-40-LED</a>	Ringlight		40 LED-HB	Segment/Intensity		Auto 100-240vac							USA	<a href="#">LED spec</a>
Techniquip	<a href="#">AN-2-80-LED</a>	Ringlight		80 LED-HB	Segment/Intensity		Auto 100-240vac							USA	<a href="#">LED spec</a>